

Ideal for temporary storage and liquid gas service, this flask line incorporates the convenience of a coated handle on the most popular wide mouth cylindrical models.

- Aluminum metal housing extends full length of Dewar for added protection
- The Dewar lip is extended to facilitate pouring
- Polyethylene stopper, specially vented, minimizes evaporation

PART#	Dewar Description	Type	Model	Protective Covering	ID mm	Inside Depth mm	OD mm	Total Height mm	Flask ml	Stopper
8900	350ml Full Base w/Metal Handle	Cylindrical, Wide Mouth	Lab Grade	Full Aluminum	70	125	86	157	–	8845
8921	1900ml Full Base w/Metal Handle	Cylindrical, Wide Mouth	Lab Grade	Full Aluminum	120	190	150	252	–	8821
8942	4300ml Full Base w/Metal Handle	Cylindrical, Wide Mouth	Lab Grade	Full Aluminum	143	305	178	365	–	8842
8945	1000ml Full Base w/Metal Handle	Cylindrical, Wide Mouth	Lab Grade	Full Aluminum	70	302	86	337	–	8845



ติดต่อสั่งซื้อได้ที่ **Pamalyne Co., Ltd.**

18/365 , . 28, ,
 , 10230
 .0-2363 8933, 0-2363 8934, 0-23963 8935 : 0-2363 8932